

UNICARE BLACK NITRILE 200 (GS900*)

Disposable examination gloves made of nitrile

Nitrile • Latex free • Powder-free • Black • Non-sterile • Single chlorinated
 Chemical protective glove EN 374 • Suitable for short-term food contact
 Protection against viruses, fungi and bacteria

Specifications:

Material: Nitrile
 Powder: Powder-free
 Ingredients: Thiuram free, no proteins contained
 Cuff: Beaded cuff
 Design: Ambidextrous
 Surface: Finger textured
 Colour: Black
 Quality assurance:
 EN 374 CAT III - Type C;
 EN ISO 21420;
 EN 455 Part 1-4;
 DIN EN ISO 13485:2021

Item Number and Size:

GS9002 | Small
 GS9003 | Medium
 GS9004 | Large
 GS9005 | Extra Large

Packaging:

200 pieces per box
 10 boxes per carton

Classification:



For more information click on the image



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical properties	EN 455-2	Median value meets the requirements
1000-ml-Watertightness test	G-1	1.5
Visual inspection for major defects	G-1	2.5
Visual inspection for minor defects	G-1	4.0

Dimensions

Item no.	GS9002	GS9003	GS9004	GS9005
Size	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	min. 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer
Palm (in the middle of the palm)	0.08 ± 0.03 mm

Physical Properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Limited, Block 1 Blanchardstown Corporate Park,
 Ballycoolin Road, Blanchardstown, D15 AKK1, Ireland
 (Notified Body 2834)

Test chemical	EN ISO 374-1:2016+A1: 2018/Typ B Permeation level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K)	6	-5.3

Hint

The above information applies to the standard version.
 Subject to changes.



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

