



SELECT Black (GT240X*)

offline washed · black · chemical tested · medical grade
medium weight · latex examination gloves

Physical Specification:

Material: Latex (Flat)
Colour: Black
Weight: 6.2 g
Thickness: 0.11 mm
Finish: Powder-free
Fit: Ambidextrous
Sterile: No

Order Codes | Sizing:

GT2401 | Extra Small
GT2402 | Small
GT2403 | Medium
GT2404 | Large
GT2405 | Extra Large

Packaging:

Each case contains 10 boxes
Each dispenser box
contains 100 pairs

Classification:



For further
information
please click
on the
image



Features

- Tested to EN374 Type B
- AQL 1.5
- Matt-black finish
- Free from powder
- Ambidextrous
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

Benefits

- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Matt black finish masks blood and ink, and reduces hand and eye strain when working on complex designs
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves"

Sectors

- Hair and beauty
- Piercing
- Tattoo



UNIGLOVES USA

UG Healthcare USA Inc.
Costa Mesa, CA 92626

(800) 328-7981
sales@ughealthcare.com
us.unigloves.com



**SAME PRODUCT
NEW LOOK**



Issue No. 1

Prepared at Gillingham
on 17/09/26