

Contact® (150*)

offline washed · chemical tested · white · medical grade
medium weight · latex · examination gloves

Physical Specification:

Material: Latex (Flat)
Colour: White
Weight: 6.2 g
Thickness: 0.11 mm
Finish: Powder-free
Fit: Ambidextrous
Sterile: No



Order Codes | Sizing:



1501 | Extra Small
1502 | Small
1503 | Medium
1504 | Large

Packaging:


Each case contains 10 boxes
Each dispenser box contains 100 pairs

Classification:

AQL 1.5 (Pinholes) EN455 (Medical Grade)




- Features**
- Washed offline
 - Tested to EN374 Type B
 - AQL 1.5
 - Free from powder
 - Tested to EN1186
 - Beaded cuff
 - Manufactured with textured palms and fingertips
 - Medium-weight construction
 - Touchscreen compatible

- Benefits**
- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
 - Resistant to a range of chemicals
 - Meets medical standards for pinholes
 - Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
 - Suitable for contact with all food types
 - Beaded cuff provides extra strength during donning and helps prevent liquid roll back
 - Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
 - Medium-weight construction provides a balance of strength and sensitivity
 - Suitable for use with mobile and tablet displays without removing gloves
 - Extremely grippy and reliable on the outside, the glove feels soft against the skin on the inside and slides effortlessly onto the hand

- Sectors**
- Surgery
 - Laboratories
 - Clinics
 - Chemical manufacture and use
 - Dentistry
 - Dental hygiene
 - Emergency services



For further information please click on the image 



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green,
Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
uk.unigloves.com

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf,
Germany
+49 (0) 2241-9323-0
unigloves@unigloves.de
de.unigloves.com

