



# LANO-E-GEL (190X\*)

moisturising · offline washed · natural · medical grade  
 medium weight · latex · examination gloves

## Physical Specifications:

Material: Premium grade  
 Latex (offline chlorinated)  
 Colour: Natural  
 Length: 240 mm  
 Weight: 6.0 g  
 Thickness: 0.11 mm  
 Finish: Powder free,  
 chlorinated  
 Fit: Ambidextrous  
 Sterile: No

## Order Codes and Size Information:

1901 | Extra small  
 1902 | Small  
 1903 | Medium  
 1904 | Large  
 1905 | Extra large

## Packaging:

Each dispenser box  
 contains 100 gloves  
 Each case contains 10  
 dispenser boxes

## Classification:



For more about  
 the **Premium**  
 range, please  
 click  
 on the  
 image



## Features


- Washed offline
- Microfilm inner coating contains a combination of natural moisturisers
- Tested to EN374
- AQL 1.5 and EN455 parts 1-4
- Free from powder
- Tested to EN1186
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

## Benefits

- Offline washing ensures ease of donning, increased comfort in use and reduced risk of allergic reactions
- Microfilm inner coating helps protect and nourish skin, and aids donning
- Resistant to a range of chemicals
- Meets medical standards including pinholes
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with all food types
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves

## Sectors

- Dentistry • Dental hygiene
- Floristry and gardening
- Hair and beauty
- Healthcare



**UNIGLOVES UK**  
 Registered in England No: 4010200  
 VAT registration 791817787

3 Ambley Green,  
 Gillingham Business Park,  
 Gillingham, ME8 0NJ, UK  
 +44 (0)800 049 6602  
 enquiries@unigloves.co.uk  
 uk.unigloves.com

**UNIGLOVES GERMANY**  
 Arzt- & Klinikbedarf  
 Handelsgesellschaft mbH

Camp-Spich-Straße 71  
 D-53842 Troisdorf,  
 Germany  
 +49 (0) 2241-9323-0  
 unigloves@unigloves.de  
 de.unigloves.com



Issue No. 1

Prepared at Gillingham  
 on 29/09/25