

Product specification sheet

Revision dated July 2026

1.	Product code	GT240*									
2.	Product description	SELECT Black									
3.	Material	Latex									
4.	Inspection parameters	Criteria		Inspection level	AQL	Standard					
		Dimensions		S-2	4.0	EN455-2					
		Physical properties		S-2	4.0	EN455-2					
		Water tight test (1,000ml)		G-1	1.5	EN455-1					
		Visual defects		G-1	2.5	In-house					
5.	Glove dimensions	Dimensions		Size (mm)							
				XS	S	M	L	XL			
		Product code	GT0021	GT0022	GT0023	GT0024	GT0025				
		Palm width	≤80	80 ± 10	95 ± 10	110 ± 10	≥ 110				
		Length	≥ 240	≥ 240	≥ 240	≥ 240	≥ 240				
6.	Glove thickness	Point of thickness measurement			Nominal single wall (mm)						
		Cuff (at beaded end)			0.08						
		Palm (at centre of palm)			0.11						
		Fingertip (at fingertip)			0.14						
7.	Physical properties	Parameters		Before aging		After aging					
		Median force at break (N)		Min 6.0		Min 6.0					
		Tensile strength (MPa)		Min 18.0		Min 14.0					
		Ultimate elongation (%)		Min 650		Min 500					
8.	Proteins	< 50 µg/g extractable latex proteins									
9.	Powder residue	Maximum 2.0 mg / glove									
10.	Cuff type	Beaded									
11.	Design	Ambidextrous									
12.	Surface	Microroughened									
13.	Colour	Black									
14.	Quality assurance	Manufactured under ISO 9001 and ISO 13485									
15.	Weight	6.1g (size medium)									
16.	Shelf life	5 years									
17.	Standards compliant with	EN455-1	EN455-2	EN455-3	EN455-4	EN ISO 21420	EN374-1	EN374-2	EN374-4	EN374-5	Food approved
		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

The above data represents the product norm, which may be changed without notice and is intended as a guide only.