



# SELECT Black (GT240X\*)

offline washed · black · chemical tested · medical grade  
medium weight · latex examination gloves

## Physical Specification:

Material: Latex (Flat)  
Colour: Black  
Weight: 6.2 g  
Thickness: 0.11 mm  
Finish: Powder-free  
Fit: Ambidextrous  
Sterile: No

## Order Codes | Sizing:


GT2401 | Extra Small  
GT2402 | Small  
GT2403 | Medium  
GT2404 | Large  
GT2405 | Extra Large

## Packaging:

Each case contains 10 boxes  
Each dispenser box contains 100 pairs

## Classification:



For further information please click on the image 



## Features

- Tested to EN374 Type B
- AQL 1.5
- Matt-black finish
- Free from powder
- Ambidextrous
- Beaded cuff
- Manufactured with textured palms and fingertips
- Medium-weight construction
- Touchscreen compatible

## Benefits

- Resistant to a range of chemicals
- Meets medical standards for pinholes
- Matt black finish masks blood and ink, and reduces hand and eye strain when working on complex designs
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Ambidextrous fit means you don't have to split a pair to replace a damaged glove
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingertips provide improved grip, particularly when working in wet or oily conditions
- Medium-weight construction provides a balance of strength and sensitivity
- Suitable for use with mobile and tablet displays without removing gloves"

## Sectors

- Hair and beauty
- Piercing
- Tattoo



**UNIGLOVES UK**  
Registered in England No: 4010200  
VAT registration 791817787

3 Ambley Green,  
Gillingham Business Park,  
Gillingham, ME8 0NJ, UK  
+44 (0)800 049 6602  
enquiries@unigloves.co.uk  
uk.unigloves.com

**UNIGLOVES GERMANY**  
Arzt- & Klinikbedarf  
Handelsgesellschaft mbH

Camp-Spich-Straße 71  
D-53842 Troisdorf,  
Germany  
+49 (0) 2241-9323-0  
unigloves@unigloves.de  
de.unigloves.com

