

UNIGLOVES® Product specification sheet Revision dated May 2022

| 1. | Product code | uct code UCHG300*G | | | | | | | | | | | | |
|-----|-----------------------------|---|-----------|-------------|-------------|---------------------|------|-----------|-------------|--------------------------|---------|------------|------------------|--|
| 2. | Product description | Unicare Household Green | | | | | | | | | | | | |
| 3. | Material | Latex | | | | | | | | | | | | |
| 4. | Inspection | Ci | riteria | | Ir | Inspection level AC | | | AQL | L Standard | | | | |
| | parameters | Dimensions | | | S-2 | | 4.0 | | | EN420 | | 120 | | |
| | | Physical properties Water tight test (1,000ml) | | | | S-2 | | | 4.0 | 4.0 | | EN455-2 | | |
| | | | | | | G-1 1 | | | 1.5 | .5 | | EN455-1 | | |
| | İ | Visual defects | | | | G-1 | | | 2.5 | 2.5 | | In-house | | |
| 5. | Glove dimensions | Size (mm) | | | | | | | | | | | | |
| | | Dimensions | | XS | | S | | M | | L | | XL | | |
| | | Product code | duct code | | ı | ICHG3002G | UC | UCHG3003G | | UCHG3004G | | UC | HG3005G | |
| | | Palm width | | - | | 95 ± 5 | | 105 ± 5 | | 110 ± 5 | 110 ± 5 | | 120 ± 5 | |
| | | Length | | - | | ≥ 300 | ≥ 3 | | 0 | ≥ 300 | | | ≥ 300 | |
| 6. | Glove thickness | Point of thickness measurement | | | | | | | | Nominal single wall (mm) | | | | |
| | | Cuff (at beaded end) | | | | | | | | 0.15 | | | | |
| | | Palm (at centre of palm) | | | | | | | | | 0. | 0.25 | | |
| | | Fingertip (at fingertip) 0.28 | | | | | | | | | | | | |
| 7. | Physical properties | Parameters | | | | Before aging | | | | After aging | | | | |
| | | Median force at break (N) | | | | N/A | | | | N/A | | | | |
| | | Tensile strength (MPa) | | | | Min 17.5 | | | | Min 13.5 | | | | |
| | | Ultimate elong | gation (| %) | | Min 600 | | | | Min 500 | | | | |
| 8. | Proteins | < 50 µg/g extractable latex proteins | | | | | | | | | | | | |
| 9. | Powder residue | n/a | | | | | | | | | | | | |
| 10. | Cuff type | Beaded | | | | | | | | | | | | |
| 11. | Design | Hand specific | | | | | | | | | | | | |
| 12. | Surface | Raised grip pattern | | | | | | | | | | | | |
| 13. | Colour | Green | | | | | | | | | | | | |
| 14. | Quality assurance | Manufactured under ISO 9001 | | | | | | | | | | | | |
| 15. | Weight | 20.0g (size medium) | | | | | | | | | | | | |
| 16. | Shelf life | | | | | | | | | | | | | |
| 17. | Standards compliant with | | 455- 2 | EN455- 3 | EN455- 4 | EN420 | EN37 | | EN374- 2 | EN374- 4 | EN | 1374- 5 | Food approved | |
| | | | | | | ✓ | | | | | | | ✓ | |