

Soft Nitril Blue (80*0)

Disposable examination gloves made of nitrile

Nitrile • Latex free • Powder-free • Blue • Non-sterile •

Suitable for short-term food contact • Chemical protective glove EN 374 •

Protection against viruses, fungi and bacteria

Specifications:

Material: Nitrile

Powder: Powder-free

Ingredients: Thiuram free,
no proteins contained

Cuff: Beaded cuff

Design: Ambidextrous

Surface: Micro-roughened

Colour: Blue

Quality assurance:

EN 455 Part 1-4;

EN 374 CAT III - Typ C;

EN ISO 21420;

DIN EN ISO 13485:2021

Item number and size:

8020 | Small

8030 | Medium

8040 | Large

8050 | Extra Large

Packaging:

100 pieces per box

10 boxes per carton

Classification:



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical properties	EN 455-2	Median value meets the requirements
1000-ml-Watertightness test	G-1	1,5
Visual inspection for major errors	G-1	2,5
Visual inspection for major errors	G-1	4,0

Dimensions

Item no.	8020	8030	8040	8050
Size	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer (mm)
Cuff (at a distance of 40 mm to the beaded edge)	min. 0.05
Palm (in the middle of the palm)	min. 0.08
Fingertips (at a distance of 15mm to the tip of the middle finger)	min. 0.09

Physical properties

Parameter	Before aging	After aging
Median force at break (N)	min. 6.0	min. 6.0

Chemical protection

The product complies with the PPE regulation (EU) 2016/425

Type examination (Module B) & Module C2

CCQS Certification Services Ltd, Blanchardstown Corporate Park, D15 AKK1, Ireland

(Notified Body 2834)

Test chemical	EN ISO 374-1:2016 Permeation level	EN ISO 374-4:2019 Degradation %
Sodium hydroxide 40% (K)	6	-5.3

Hint

The above information applies to the standard version.

Subject to changes.



UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
 D-53842 Troisdorf
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de

SUBSIDIARY OF THE

UNIGLOVES®
 UG HEALTHCARE GROUP
 ughealthcarecorporation.com