

Disposable paper Napkins (51** *)

Absorbent and liquid-impermeable napkin

Cellulose wadding • in different colours

Specifications:

Material: Cellulose wadding

Colour and Item-No.:

burgundy	5102
cedro	5103
magenta	5107
mint green	5120
light blue	5121
white	5122
yellow	5123
pink	5124
black	5125
orange	5126
dark blue	5128
lilac	5129

Packaging:

50 pieces per bag
10 bags per carton

For more
information
click on
the image



Dimension

Art.-No.	colour	Length	Width
5102	burgundy	33 cm	45 cm
5103	cedro	33 cm	45 cm
5107	magenta	33 cm	45 cm
5120	mint green	33 cm	45 cm
5121	light blue	33 cm	45 cm
5122	white	33 cm	45 cm
5123	yellow	33 cm	45 cm
5124	pink	33 cm	45 cm
5125	black	33 cm	45 cm
5126	orange	33 cm	45 cm
5128	dark blue	33 cm	45 cm
5129	lilac	33 cm	45 cm



Thickness

Cellulose wadding, 16.5 g/m² per ply, polyethylene 9µm +/- 10%

Sterilization & Durability

Disposable product, not suitable for sterilization, 5 years shelf life.

Storage

Store in a dry and clean place, away from heat sources.

Hint

The above information applies to the standard version.
Subject to changes.