

Expert Plus (12**)

Disposable Latex Examination Gloves

Latex • Powder free • White • Sterile • Anatomical fit

Specifications:

Material: Natural High Grade Rubber Latex
 Powder: Powder free (< 2 mg/glove)
 Protein content: < 50 µg/g extractable latex proteins
 Cuff: beaded cuff
 Design: Hand-specific
 Colour: White
 Quality control:
 EN 455 Teil 1-4;
 DIN EN ISO 13485:2021

Item number and size:

1260 | 6.0
 1265 | 6.5
 1270 | 7.0
 1275 | 7.5
 1280 | 8.0
 1285 | 8.5

Packaging:

50 pair per dispenser
 4 dispensers per carton

Classification:



Test parameter

Test criteria	Level	AQL
Dimensions	EN 455-2	Median value meets the requirements
Physical properties	EN 455-2	Median value meets the requirements
1000-ml-Watertightness test	G-1	0.65
Visual inspection for major errors	G-1	0.65
Visual inspection for minor errors	G-1	2.5

Dimensions

Item-no.	1260	1265	1270	1275	1280	1285
Size	6.0	6.5	7.0	7.5	8.0	8.5
Width	77 ± 5 mm	83 ± 5 mm	89 ± 5 mm	95 ± 5 mm	102 ± 6 mm	108 ± 6 mm
Length	≥ 295 mm	≥ 295 mm	≥ 295 mm	≥ 295 mm	≥ 295 mm	≥ 295 mm

Material thickness

Measuring	One layer (mm)
Cuff (at a distance of 25 mm to the beaded edge)	0.18 ± 0.02
Palm (in the middle of the palm)	0.22 ± 0.02
Fingertips (at a distance of 15mm to the tip of the middle finger)	0.23 ± 0.02

Physical properties

Parameter	Before aging	After aging
Force at Break (N)	min. 9.0	min. 9.0

Quality assurance

Comply the requirements of ASTM D3577-91, EN 455 part 1-4
 Manufactured according GMP and EN ISO 9001 and EN 13485 quality assurance system

Hint

The above information apply to the standard version.
 Subject to change.



TERANG NUSA Sdn. Bhd.
 1 Jalan 8
 Pengkalan Chepa 2 Industrial Zone
 16100, Kota Bharu, MALAYSIA



ULMA INTERNATIONAL GmbH
 Pfaffenweg 35
 89231 Neu-Ulm, Germany
 MADE IN MALAYSIA



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green,
 Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 uk.unigloves.com

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf,
 Germany
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 de.unigloves.com



Issue No. 1.2

Prepared at Gillingham
 on 07/08/2026