

PROFIL PLUS SMALL (S2005 X)

Surgical Face Mask with elastic band

Soft, round, non-irritating elastic bands •
 3-layer non-fiber-glass filter medium • latex free •
 integrated nickel-free nose bridge

Item-No. and colour:

S2005 blue	Blue
S2005 black	Black
S2005 pink	Pink
S2005 white	White

Packaging:

50 Pieces per Box
 10 Boxes per Carton



Xiamen EAN Protecting
 Products Co., Ltd.
 No. 79-85 ZhenNan Rd3,
 TongAn District, Xiamen,
 Fujian, 36110, China
 Made in China



Unigloves (Germany)
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH
 Camp-Spich-Straße 71
 D-53842 Troisdorf
 Germany

Classification:



For more
 Information
 click on
 the
 Image



Intended use according

The product is non-sterile mask, which is used in the patient's environment in order to protect the patient and the user from contamination.

Bacterial filtration efficiency (BFE)

BFE $\geq 98\%$ according to EN 14683, Type IIR.

Dimensions

of ready for use mask

Width	95 ± 2 mm
Length	145 ± 2 mm

Outer layer

Material	Thickness
100 % Polypropylene non-woven	15 g/m^2

Middle layer

Material	Thickness
„Melt Blown Filter“	25 g/m^2

Inner layer

Material	Thickness
100 % Polypropylene non-woven	30 g/m^2

Elastic Band

Material	Length	Width
braided latex-free elastic cord	175 ± 5 mm	$3,5 \pm 0,2$ mm

Nose bridge

Material	Length	Width
wire coated with plastic	$110 \text{ mm} \pm 2$ mm	$2,5 \text{ mm} \pm 0,5$ mm

Quality assurance

All materials used are free of irritants and odorless, do not contain glass fibers and are suitable for medical and dental purposes.

Hint

The above information applies to the standard version. Subject to changes.



UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

📍 Camp-Spich-Straße 71
 D-53842 Troisdorf
 ☎ +49 (0) 2241-9323-0
 ✉ unigloves@unigloves.de
 🌐 unigloves.de

