



Azo™ Alcohol Efficacy Summary

Micro-organism	European Standard	Contact Time
Bactericidal		
Acinetobacter baumannii (MDR)	EN 13727	15 seconds
Campylobacter jejuni	EN 1276	15 seconds
Enterococcus hirae	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
Escherichia coli	EN 1276	15 seconds
	EN 13697	60 seconds
Klebsiella pneumoniae	EN 13727	15 seconds
Legionella pneumophila	EN 1276	15 seconds
Listeria monocytogenes	EN 1276	15 seconds
Pseudomonas aeruginosa	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 13697	60 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
Salmonella enterica	EN 1276	15 seconds
Staphylococcus aureus	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 13697	60 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
MRSA	EN 13727	15 seconds
Vancomycin resistant enterococcus (VRE)	EN 13727	15 seconds
Mycobacterial		
Mycobacterium terrae	EN 14348	60 seconds

Micro-organism	European Standard	Contact Time
Fungicidal and Yeastocidal		
Candida albicans	EN 13624	15 seconds
	EN 1650	15 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
Aspergillus brasiliensis	EN 14562	60 seconds
Virucidal		
Bovine viral diarrhoea virus (Hep C)	EN 14476	15 seconds
Influenza A H1N1	EN 14476	15 seconds
Vaccinia virus	EN 14476	15 seconds
Feline coronavirus	EN 14476	15 seconds
Murine Norovirus	EN 14476	60 seconds
Duck Hepatitis Virus (Hep B)	ASTM E1052	15 seconds

Azo™ 70% IPA wipes and sprays are compliant to biocidal product regulation 528/2012 and are manufactured to ISO 13485 standards. **Use biocides safely. Always read the label and product information before use.**

For use in care environments, including:



Hospital Wards



Laboratories



Doctors Clinics



Cleanrooms



Operating Theatres



Aesthetics Clinics

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
 +44 (0)1204 529494 hello@vernacare.com vernacare.com

Before use, please familiarise yourself with the Instructions for Use and ensure you follow any precautions and warning stated on the packaging. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
 © Copyright Vernacare International Limited 2025. All rights reserved. IC203 REV01

