

TECHNICAL DATA SHEET

128AA080 Conti AQUAPURE Washbowl

CONTI®

AQUAPURE WASHBOWL

Name: Conti AQUAPURE Washbowl
Product Code: 128AA080
Quantity: 80 units per carton



General Information

Capacity: 4.0 litres max (2.0 litres recommended)
Material: Recycled cellulose fibres
Weight: 51g (+/- 5g)
Colour: Grey
Unit Dimensions: 341x259x93mm
Case Dimensions: 570x340x280mm
Expiry: 1 Years from date of manufacture

Compliance & Performance

- Tested in accordance with PAS 29:1999
- ISO 9001:2008 & ISO 14001:2014 Certified
- Guaranteed fluid & detergent retention for 2 hours.
- Resistance temperature: 45° C

GTIN Code: 5 010058 000698



Key Features

- Fluorocarbon (PFAS) Forever Chemical Free Product*
- Made for patient bathing and cleaning
- Designed to hold detergent with warm water
- Fully disposable, for use with Vernacare macerators

*No PFAS has been added to the product during manufacture, though trace amounts may be present as a result of environmental distribution.



Protecting Safety & Dignity. Pioneering Sustainable Manufacturing.

/ infection prevention / clinical waste management / surgical solutions

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

Vernacare International Limited. Registered in England and Wales. Company No: 06449040.
Registered office: 1 Western Avenue, Buckshaw Village, Chorley, England, PR7 7NB.

Vernacare