



# Azo™ Universal Efficacy Summary



| Micro-organism                                     | European Standard | Contact Time |
|----------------------------------------------------|-------------------|--------------|
| <b>Bactericidal</b>                                |                   |              |
| Acinetobacter baumannii (MDR)                      | EN 13727          | 30 seconds   |
| Enterococcus faecalis                              | EN 13727          | 30 seconds   |
|                                                    | EN 13697          | 60 seconds   |
| Enterococcus hirae                                 | EN 13727          | 30 seconds   |
|                                                    | EN 1276           | 60 seconds   |
|                                                    | EN 16615          | 2 minutes    |
| Escherichia coli                                   | EN 13727          | 30 seconds   |
|                                                    | EN 1376           | 60 seconds   |
|                                                    | EN 13697          | 60 seconds   |
| Klebsiella pneumoniae                              | EN 13727          | 30 seconds   |
|                                                    | EN 1276           | 60 seconds   |
| Listeria monocytogenes                             | EN 13727          | 30 seconds   |
| Pseudomonas aeruginosa                             | EN 13727          | 30 seconds   |
|                                                    | EN 1276           | 60 seconds   |
|                                                    | EN 13697          | 60 seconds   |
|                                                    | EN 16615          | 2 minutes    |
| Salmonella enterica ser. Typhimurium               | EN 1276           | 60 seconds   |
| Staphylococcus aureus                              | EN 13727          | 30 seconds   |
|                                                    | EN 1276           | 60 seconds   |
|                                                    | EN 16615          | 2 minutes    |
| Methicillin-resistant Staphylococcus aureus (MRSA) | EN 13727          | 30 seconds   |
|                                                    | EN 1276           | 60 seconds   |
| <b>Mycobacterial</b>                               |                   |              |
| Mycobacterium terrae                               | EN 14348          | 5 minutes    |
| <b>Fungicidal and Yeastcidal</b>                   |                   |              |
| Candida albicans                                   | EN 16615          | 2 minutes    |
|                                                    | EN 13624          | 5 minutes    |
| Candida auris                                      | EN 13624          | 2 minutes    |
| <b>Virucidal</b>                                   |                   |              |
| Human coronavirus                                  | EN 14476          | 30 seconds   |
| Vaccinia virus                                     | EN 14476          | 30 seconds   |
| Feline coronavirus                                 | EN 14476          | 60 seconds   |
| Adenovirus type 5                                  | EN 14476          | 5 minutes    |
| Murine Norovirus                                   | EN 14476          | 5 minutes    |

**Medical Device:**  
For professional use only  
within healthcare environments,  
including:



**Care  
Homes**



**Hospital  
Wards**



**Doctors  
Clinics**



**Dental  
Practices**



**Pharmacies**



**Podiatry  
Clinics**

Class IIa medical device compliant – Azo™ Universal wipes are compliant to medical device directive 93/42/EEC, and are manufactured to ISO 13485 standards.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK  
+44 (0)1204 529494 hello@vernacare.com vernacare.com

Before use, please familiarise yourself with the Instructions for Use and ensure you follow any precautions and warning stated on the packaging. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.  
© Copyright Vernacare International Limited 2025. All rights reserved. IC202 REV01



**Vernacare®**