

The azo logo consists of the word "azo" in a white, lowercase, sans-serif font, with a registered trademark symbol (®) to its upper right. It is set against a dark blue, rounded, pill-shaped background.

Universal Plus

Designed to deliver more with
The Power of Plus



What Does Plus Mean?

Plus⁺ Protection

Protecting what matters most with dual action cleaning and disinfection.

Plus⁺ Confidence

Confidence in every clean with proven efficacy.

Plus⁺ Convenience

More time for care and less time preparing.

Plus⁺ Flexibility

A choice for different clinical environments with performance in the format you need.

The Vernacare logo features a stylized blue wave icon above the word "Vernacare" in a dark blue, sans-serif font, with a registered trademark symbol (®) to its upper right.

Fast, Proven Efficacy

Full efficacy data available on request.

Organism	EN Standard	Contact Time
Bactericidal	EN 16615	2 minutes
	EN 13727	30 seconds
Limited Spectrum Virucidal	EN 14476	1 minute
Yeasticidal	EN 13624	30 seconds
	EN 16615	2 minutes
Tuberculocidal	EN 14348	2 minutes

Ordering Information

Product Code	Description	Pack Size	Case Breakdown
AUW50FW	Azo® Universal Plus Dual Action – Cleaning and Disinfection Wipes	50	50 wipes per pack x 24 packs per case
AUW200FW	Azo® Universal Plus Dual Action – Cleaning and Disinfection Wipes	200	200 wipes per pack x 6 packs per case



To find out more, speak to your Commercial Account Manager, or **visit www.vernacare.com** today.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

Before use, please familiarise yourself with the Instructions for Use and ensure you follow any precautions and warning stated on packaging. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
© Copyright Vernacare International Limited 2026. All rights reserved. IC173 REV00



Vernacare®