



Azo™ Alcohol Efficacy Summary

Micro-organism	European Standard	Contact Time
Bactericidal		
Acinetobacter baumannii (MDR)	EN13727	15 seconds
Campylobacter jejuni	EN 1276	15 seconds
Enterococcus hirae	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
Escherichia coli	EN 1276	15 seconds
	EN 13697	60 seconds
Klebsiella pneumoniae	EN13727	15 seconds
Legionella pneumophila	EN 1276	15 seconds
Listeria monocytogenes	EN 1276	15 seconds
Pseudomonas aeruginosa	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 13697	60 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
Salmonella enterica	EN 1276	15 seconds
Staphylococcus aureus	EN 13727	15 seconds
	EN 1276	15 seconds
	EN 13697	60 seconds
	EN 14561	60 seconds
	EN 16615	60 seconds
MRSA	EN 13727	15 seconds
Vancomycin resistant enterococcus (VRE)	EN 13727	15 seconds
Mycobacterial		
Mycobacterium terrae	EN 14348	60 seconds

Micro-organism	European Standard	Contact Time
Fungicidal and Yeastocidal		
Candida albicans	EN 13624	15 seconds
	EN 14561	15 seconds
	EN 16615	60 seconds
	EN 1650	60 seconds
Aspergillus brasiliensis	EN 14562	60 seconds
Virucidal		
Bovine viral diarrhoea virus (Hep C)	EN 14476	15 seconds
Influenza A H1N1	EN 14476	15 seconds
Vaccina virus	EN14476	15 seconds
Feline coronavirus	EN14476	15 seconds
Murine Norovirus	EN14476	60 seconds
Duck Hepatitis Virus (Hep B)	ASTM E1052	15 seconds

Class IIa medical device compliant – Azo™ Alcohol Equipment wipes are compliant to medical device directive 93/42/EEC, and are manufactured to ISO 13485 standards.

Medical Device - For professional use in care environments, including:



Hospital Wards



Doctors Clinics



Dental Practices



Podiatry Clinics

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

Before use, please familiarise yourself with the Instructions for Use and ensure you follow any precautions and warning stated on the packaging. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
© Copyright Vernacare International Limited 2025. All rights reserved. IC204 REV01



Vernacare®