

Clinical Waste Management Range Card

NHS – England

Clinical Waste Management



Vernacare®

Contents

Sharpsafe®

- 4 Pocket Range 0.2L / 0.45L / 0.6L
- 5 Small Range 1L / 1.8L / 2L / 3L
- 6 Medium Range 5L / 7L
- 7 Large Range 9L / 13L / 24L / 30L
- 8 Specialist 7.5L / H9L / 25XL

Clinisafe®

- 9 Plastic Range 3L / 8.5L / 11.5L / 30L / 50L / 60L / UN3549 60L
- 11 Cardboard Range 5L / 12L / 25L / 50L

12 **Accessories**

Brackets / Stands / Trays / Trolleys

Welcome to Vernacare Clinical Waste Management

Pioneering a safe and sustainable Clinical Waste Management solution.

Here at Vernacare, we are passionate about clinical waste management and understand the healthcare operations needed to carefully handle, segregate and properly dispose of clinical waste.

Over 40 years ago, we pioneered the first purpose-built sharps container and we continue to innovate. We've introduced the world's first full range of sustainable, single-use clinical waste containers and established the Near Patient Sharps Disposal system for safe disposal of sharps, combating needlestick injuries.

Our clinical waste management portfolio offers a comprehensive range of products to dispose of sharps, pathological waste, pharmaceuticals, and more.

We are sustainable but never at the compromise of safety.

We're committed to providing safe, sustainable products that not only protect patients and healthcare professionals, but also reduce our impact on the environment.

We are committed to helping you reduce your carbon footprint. Our Sharpsafe® containers have a 94% reduction in their carbon footprint when compared to containers using virgin plastic. Our range combines high-performance sharps waste management with our eco-conscious manufacturing. By switching our manufacturing over to low carbon recycled materials, we are helping to support healthcare facilities to achieve their carbon Net Zero targets.



Up to 100% recycled material in the base of the containers



94% reduction in carbon footprint compared to virgin containers



Zero waste to landfill utilising high-grade reprocessed material, efficient use of energy and natural resources



Reducing our CO2 emissions through solar panels at our manufacturing site



Pocket Range 0.2L / 0.45L / 0.6L

The pocket range is designed to fit neatly into a bag or pocket for safe sharps disposal whilst in the community.

Key Features:

- ✓ Recycled Material
- ✓ Temporary and Permanent Closure Feature
- ✓ Compact Size

Imagery not to scale.



Sharpsafe® Pocket Range				Yellow		Orange		Purple	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code
Pocket	0.2L	100	Petal Opening	51731430	FNC85469	-	-	-	-
Pocket	0.2L	100	Clear Opening	51732430	FNC85480	-	-	-	-
Pocket	0.2L	100	Pen Opening	51733430	FNC85471	-	-	-	-
Pocket	0.45L	100	Petal Opening	51701430	FSL2304	-	-	-	-
Pocket	0.45L	100	Clear Opening	51702430	FNC85075	-	-	-	-
Pocket	0.45L	100	Pen Opening	51703430	FNC85076	-	-	-	-
Pocket	0.6L	100	Clear Opening	51502430	FSL2300	51502410	FSL2352	51502420	FSL2319

Accessories

Description	Qty	Code	NPC Code
0.6L Stand	10	41510000	-
0.45/0.6L Wall Bracket	10	42500000	-

Small Range 1L / 1.8L / 2L / 3L

The small range is perfectly designed to help with the safe disposal of sharps in the community and acute settings. Compatible with a range of accessories to help minimise the risk of sharps injuries at the point of use.

Key Features:

- ✓ Recycled Material
- ✓ Temporary and Permanent Closure Feature
- ✓ In-mould Label

Imagery not to scale.



Sharpsafe® Small Range				Yellow		Orange		Purple	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code
Small	1L	100	-	51602430	FSL2347	51602410	FSL2325	51602420	FSL2320
Small	1.8L	75	-	50502430	FNC85063	-	-	50502420	FNC85064
Small	2L	50	Safety Shield	51405430	FSL2317	51405410	FSL2323	51405420	FSL2299
Small	3L	50	Safety Shield	51455430	FSL2328	51455410	FSL2339	51455420	FSL2311

Accessories

Description	Qty	Code	NPC Code
1L Stand (Printed)	52	42530000	FNC85049
1L/1.8L/2L/3L Stand (Glue Tabs)	26	42050000	FNC85070
1L/1.8L/2L/3L NPSD Tray	10	43000000	FNC85061
1L/1.8L/2L/3L Triple Tray	10	43010000	-
2L/3L Backing Board Bracket	1	42870000	-
SafeClip Wall Bracket	10	45000000	FNC85041
1.8L Wall Safe	1	42000005	-

sharp safe® With DUAL PROTECT PRO™

Medium Range 5L / 7L

The new medium range is designed to maximise capacity in a small footprint for acute areas. With innovative safety features such as the Dual Protect Pro™ these containers will help minimise the risk of sharps injuries, whilst lowering their impact on the environment.

Key Features:

- ✓ 100% Recycled Plastic Base
- ✓ Integrated Safety Shield
- ✓ Increased Capacity

Imagery not to scale.



Sharpsafe® Medium Range				Yellow		Orange		Purple		Blue	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Mid	5L	48	Dual Protect Pro™	51025430	FNC85331	51025410	FNC85327	51025420	FNC85334	51025440	FNC85333
Mid	7L	48	Dual Protect Pro™	51125430	FNC85332	51125410	FNC85330	51125420	FNC85329	51125440	FNC85328

Accessories

Description	Qty	Code	NPC Code
SafeClip Wall Bracket	10	45000000	FNC85041
5L/7L Stand	10	42570000	FNC85477
5L/7L Backing Board Bracket	1	42720000	FNC85492

Specially designed for bespoke nuclear medicine radiation shields, these containers arrive preassembled

Range	Size	Qty	Description	Code	NPC Code
Mid	4L	12	Safety Shield	51005431	FNC85494
Mid	7L	8	Safety Shield	51105431	FNC85493

Large Range 9L / 13L / 24L / 30L

The large range containers each have a wide aperture for the disposal of larger sharps items. This range offers a protected access feature which prevents hands and unwanted items from entering the aperture.

Key Features:

- ✓ Recycled Material
- ✓ Protected Access Available
- ✓ In-mould Label

Imagery not to scale.



Sharpsafe® Large Range				Yellow		Orange		Purple		Blue	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Large	9L	20	-	51172430	FSL2324	51172410	FSL2342	51172420	FSL2356	-	-
Large	9L	20	Protected Access	51171430	FSL2306	51171410	FNC85046	51171420	FNC85060	-	-
Large	13L	20	-	51152430	FSL2349	51152410	FSL2338	51152420	FSL2322	51152440	FNC85470
Large	13L	20	Protected Access	51151430	FSL2344	51151410	FSL2334	51151420	FSL2330	-	-
Large	24L	15	-	51202431	FSL2327	51202411	FSL2353	51202421	FSL2333	51202441	FNC85085
Large	24L	10	Protected Access	51201430	FSL2301	51201410	FNC85066	51201420	FSL2329	-	-
Large	30L	15	-	51802431	FSL2321	51802411	FSL2345	51802421	FSL2346	51802441	FSL2332
Large	30L	10	Protected Access	51801430	FSL2305	51801410	FNC85056	51801420	FSL2308	51801440	-

Accessories

Description	Qty	Code	NPC Code
9L/13L Wall Loop	1	42150000	FNC85051
24L/30L Wall Loop	1	42200000	FNC85083
9L/13L Backing Plate Bracket	1	42410000	-
24L/30L Backing Plate Bracket	1	42210000	FNC85040
9L/13L Stand	20	42070000	FNC85087
9L/13L/24L/30L Procedure Trolley	1	44000000	FNC85081
9L/13L/24L/30L Pedestal Trolley	1	44800000	FNC85054



Specialist 7.5L / H9L / 25XL

The specialist range is designed to store and dispose of larger sharp instruments within a theatre environment. For guarantee of stability all containers in this range must be used with their accompanying reusable accessory.

Key Features:

- ✓ Recycled Material
- ✓ Temporary and Permanent Closure
- ✓ Suitable for Theatre

Imagery not to scale.



Sharpsafe® Specialist Range				Yellow		Orange	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code
Specialist	7.5L	20	-	51902430	FNC85057	-	-
Specialist	7.5L	20	Protected Access	51901430	FNC85080	-	-
Specialist	H9L	40	Safety Shield	51115430	FSL2313	-	-
Specialist	25XL	10	-	51952430	FSL2354	51952410	FNC85062
Specialist	25XL	10	Protected Access	51951430	FSL2343	51951410	-

Accessories			
Description	Qty	Code	NPC Code
7.5L Wall Loop	1	42750000	FNC85065
High 9L Stand	10	42520004	-
25XL Trolley	1	44710000	FNC85055

Plastic Range 3L / 8.5L / 11.5L / 30L / 50L / 60L

Available in small, medium and large sizes, the Clinisafe® container range offers a solution for the disposal of clinical waste within most healthcare settings.

Key Features:

- ✓ Recycled Material
- ✓ Hermetically Sealed
- ✓ Waterproof

Imagery not to scale.



Clinisafe® Plastic Range			Yellow		Orange		Purple		Blue		Red	
Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
3L	20	-	50102430	FSL85007	-	-	-	-	-	-	50102415	FSL2348
8.5L	20	-	51912430	FNC85082	51912410	FNC85476	-	-	-	-	51912415	FNC85059
11.5L	20	-	51922430	FSL2309	51922410	FNC85475	-	-	-	-	51922415	FSL2340
30L	10	Flat Lid	51850000	FSL2314	-	-	51850420	FNC85078	51850540	FSL2316	51850415	FSL2336
30L	10	Porthole Lid	51850432	FSL2331	-	-	51850422	-	51850542	-	-	-
50L	10	Flat Lid	51840000	-	-	-	51840420	-	51840540	-	51840415	-
50L	10	Porthole Lid	51840432	-	-	-	51840422	-	51840542	-	-	-
60L	10	Flat Lid	51860000	FSL2315	51860410	FNC85473	51860420	FSL2341	51860540	FSL2298	51860415	FSL2331
60L	10	Porthole Lid	51860432	-	51860412	FNC85484	51860422	-	51860542	-	-	-

Accessories

Description	Qty	Code	NPC Code
8/11L Stand	1	42070000	FNC85087
8/11L Wall Loop	1	42150000	FNC85051
8/11L Backing Plate Bracket	1	42410000	-
30L/50L/60L Trolley	1	Coming Soon	Coming Soon
Absorbent Pads	2000	026100	MXR85000
Clinisafe Bags	2000	020197	-



UN3549 60L

Our 60L Clinisafe® UN3549 Container is fully certified to hold Category A waste during HCID outbreaks. Suitable for; storing soft waste sharps containers that have been bagged according to ADR guidance and placed into yellow bulk-approved bags prior to disposal in the container. An extra large label makes safety information easy to read and provides a permanent audit trail.

Key Features:

- ✓ Rubber Hermetic Seal
- ✓ Waterproof
- ✓ Extra Large Self Adhesive Label

Imagery not to scale.



Clinisafe® UN3549

Size	Qty	Description	Yellow	
			Code	NPC Code
60L	10	Flat Lid	51860001	Coming Soon

Accessories

Description	Qty	Code	NPC Code
30L/50L/60L Trolley	1	Coming Soon	Coming Soon
Absorbent Pads	2000	026100	MXR85000

Cardboard Range 5L / 12L / 25L / 50L

The Clinisafe® cardboard range offers a sustainable solution to the management of clinical waste.

Key Features:

- ✓ Recycled Content Cardboard
- ✓ Easy to Assemble
- ✓ Arrives Flat-packed for Easy Storage

Imagery not to scale.



Clinisafe® Cardboard Range			Yellow		Orange		Purple		Blue		Gypsum	
Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
5L	10/800	-	20051001	Coming Soon	20053001	Coming Soon	-	-	20055001	Coming Soon	-	-
12L	10/800	-	20501001	FNC85336	20503001	FNC85337	-	-	20505001	FNC85338	-	-
25L	10/800	-	21001001	FNC85339	21003001	FNC85340	-	-	21005001	FNC85341	21007001	FNC85489
50L	10/400	-	-	-	21503001	FSL2125	21502001	-	21505001	FSL2133	21507001	FNC85342
50L	10/400	Tall	22001001	FSL2128	-	-	-	-	-	-	-	-

Clinisafe® Cardboard Range			Offensive Waste		Accessories			
Size	Qty	Description	Code	NPC Code	Description	Qty	Code	NPC Code
5L	10/800	-	20056001	Coming Soon	25L Clinisafe® Cardboard Plinth	1	44100004	FNC85045
25L	10/800	-	21006001	Coming Soon	50L Short Clinisafe® Cardboard Plinth	1	44100001	FNC85079
50L	10/800	Short	21506001	Coming Soon	50L Tall Clinisafe® Cardboard Plinth	1	44100005	FNC85044


Accessories Brackets / Stands

Imagery not to scale.

Brackets			
Description	Qty	Code	NPC Code
0.6L Wall Bracket	10	42500000	-
7.5L Wall Loop	1	42750000	FNC85065
8L/9L/11L/13L Wall Loop	1	42150000	FNC85051
24L/30L Wall Loop	1	42200000	FNC85083
2L/3L Backboard Bracket	1	42870000	-
5L/7L Backboard Bracket	1	42720000	FNC85492
8L/9L/11L/13L Backboard Bracket	1	42410000	-
24L/30L Backboard Bracket	1	42210000	FNC85040
SafeClip Wall Bracket	10	45000000	FNC85041



Stands			
Description	Qty	Code	NPC Code
0.6L Stand	10	41510000	-
1L Stand	52	42530000	FNC85049
1L/1.8L/2L/3L Stand (Glue Tabs)	26	42050000	FNC85070
5L/7L Stand	10	42570000	FNC85477
8L/9L/11L/13L Stand	20	42070000	FNC85087
High 9L Stand	10	42520004	-
25L Clinisafe® Cardboard Plinth	1	44100004	FNC85045
50L Short Clinisafe® Cardboard Plinth	1	44100001	FNC85079
50L Tall Clinisafe® Cardboard Plinth	1	44100005	FNC85044



Accessories Trays / Trolleys / Additional

Trays

Description	Qty	Code	NPC Code
1L/1.8L/2L/3L NPSD Tray	10	43000000	FNC85061
1L/1.8L/2L/3L Triple Tray	10	43010000	-
1L/1.8L/2L/3L Blood Tray	12	43050000	-



Trolleys

Description	Qty	Code	NPC Code
9L/13L/24L/30L Procedure Trolley	1	44000000	FNC85081
25XL Trolley	1	44710000	FNC85055
9L/13L/24L/30L Pedestal Trolley	1	44800000	FNC85054
SafeClip Trolley Adapter – Round 25mm	5	45000003	FNC85069
SafeClip Trolley Adapter – Square 25mm	5	45000004	FNC85088
SafeClip Trolley Adapter – Medi-Rail*	5	45000005	FNC85042



Additional

Description	Qty	Code	NPC Code
6L Community Specimen Bag	1	50010000	FNC85072
19L Community Specimen Bag	1	50020000	FNC85053
Absorbent Pads C135 – 150 x 150mm	2000	026100	MCR8 5000
Clinisafe® Bags	2000	020197	-

*Additional Medi-Rail adapter sizes available. Please contact your Key Account Manager for more information.

Notes

Notes

Infection Prevention

Patient Hygiene | Conti®

Environmental Cleaning | Azo™

Human Waste Disposal | Vortex Air | Vernacare Medical Pulp

Clinical Waste Management

Sharps Containers | Sharpsafe®

Clinical Waste Containers | Clinisafe®

Surgical Solutions

Women's Health | Ultraspec® | EcoGold® | Instrapac®

Laparoscopic | bert® | hubert | nubert | albert

Woundcare | Comfi™ | Cestra® Ultra four | Skintact®

Swabs & Gauze | Progauze

ENT | Instrapac® | Exmoor®

Podiatry | Instrapac®

Dental | Instrapac®



1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. © Copyright Vernacare International Limited 2025. All rights reserved. CWM152 REV00