

# Clinical Waste Management Range Card

NHS – Scotland



WARNING: DO NOT FILL ABOVE THE LINE



# Welcome to Vernacare Clinical Waste Management

## Pioneering a safe and sustainable Clinical Waste Management solution.

Here at Vernacare, we are passionate about clinical waste management and understand the healthcare operations needed to carefully handle, segregate and properly dispose of clinical waste.

Over 40 years ago, we pioneered the first purpose-built sharps container and we continue to innovate. We've introduced the world's first full range of sustainable, single-use clinical waste containers and established the Near Patient Sharps Disposal system for safe disposal of sharps, combating needlestick injuries.

Our clinical waste management portfolio offers a comprehensive range of products to dispose of sharps, pathological waste, pharmaceuticals, and more.

## We are sustainable but never at the compromise of safety.

We're committed to providing safe, sustainable products that not only protect patients and healthcare professionals, but also reduce our impact on the environment.

We are committed to helping you reduce your carbon footprint. Our Sharpsafe® containers have a 94% reduction in their carbon footprint when compared to containers using virgin plastic. Our range combines high-performance sharps waste management with our eco-conscious manufacturing. By switching our manufacturing over to low carbon recycled materials, we are helping to support healthcare facilities to achieve their carbon Net Zero targets.

### Key stats:



Up to 100% recycled material in the base of the containers



94% reduction in carbon footprint compared to virgin containers



Zero waste to landfill utilising high-grade reprocessed material, efficient use of energy and natural resources



Reducing our CO2 emissions through solar panels at our manufacturing site





sharpsafe 2L

**DANGER**  
CONTAMINATED  
SHARPS ONLY

WARNING: DO NOT FILL ABOVE THE LINE

TO BE INCINERATED

6

Vernacare



# Pocket Range

0.6L



Imagery not to scale

The pocket range is designed to fit neatly into a bag or pocket for safe sharps disposal whilst in the community.

**Key Features:** ☒ Recycled Material ☒ Temporary and Permanent Closure Feature ☒ Compact Size

Sharpsafe® Pocket Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NHS Code	Code	NHS Code	Code	NHS Code	Code	NHS Code
Pocket	0.6L	100	Clear Opening	51502410	276462	51502420	276455	-	-	-	-



## Small Range

1L / 1.8L / 2L / 3L



Imagery not to scale

The small range is perfectly designed to help with the safe disposal of sharps in the community and acute settings. Compatible with a range of accessories to help minimise the risk of sharps injuries at the point of use.

**Key Features:** ☒ Recycled Material ☒ Temporary and Permanent Closure Feature ☒ In-mould Label

Sharpsafe® Small Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Small	1L	100		51602410	276493	51602420	276486	51602220	276813	51602280	-
Small	1.8L	75		-		50502420	-	-		-	
Small	2L	50	Safety Shield	51405410	276523	51405420	276516	-		-	
Small	3L	50	Safety Shield	51455410	276547	51455420	-	51455220	277995	51455280	-

# sharpsafe® With DUAL PROTECT **PRO**™

## Medium Range 5L / 7L



Imagery not to scale

The new medium range is designed to maximise capacity in a small footprint for acute areas. With innovative safety features such as the Dual Protect Pro™ these containers are designed to help minimise the risk of sharps injuries, whilst lowering their impact on the environment.

**Key Features:**
 ✓ 100% Recycled Plastic Base
 ✓ Star Disconnection Feature
 ✓ Increased Capacity

Sharpsafe® Medium Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Medium	5L	48	Dual Protect Pro™	51025410	-	51025420	-	51025220	-	51025280	-
Medium	7L	48	Dual Protect Pro™	51125410	285983	51125420	285990	51125220	285969	51125280	285976

## Large Range

9L / 13L / 24L / 30L



Imagery not to scale

The large range containers each have a wide aperture for the disposal of larger sharps items. This range offers a protected access feature which prevents hands and unwanted items from entering the aperture.

**Key Features:** ☒ Recycled Material ☒ Protected Access Available ☒ In-mould Label

Sharpsafe® Large Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Large	9L	20		51172410	276769	51172420	-	-	-	-	-
Large	9L	20	Protected Access	51171410	-	51171420	-	-	-	-	-
Large	13L	20		51152410	276677	51152420	276660	51152220	276844	51152280	279507
Large	13L	20	Protected Access	51151410	276622	51151420	-	-	-	-	-
Large	24L	15		51202411	276691	51202421	276684	51202221	284382	-	-
Large	24L	10	Protected Access	51201410	-	51201420	-	-	-	-	-
Large	30L	15		51802411	276738	51802421	276721	51802221	276851	-	-
Large	30L	10	Protected Access	51801410	276707	51801420	-	-	-	-	-



sharpsafe®

Specialist Range

25XL



Imagery not to scale

The specialist range is designed to store and dispose of larger sharp instruments within a theatre environment. For guarantee of stability all containers in this range must be used with their accompanying reusable accessory.

**Key Features:** ☒ Recycled Material ☒ Temporary and Permanent Closure ☒ Suitable for Theatre

Sharpsafe® Specialist Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Specialist	25XL	10		51952410	-	-	-	-	-	51952280	-
Specialist	25XL	10	Protected Access	51951410	-	-	-	-	-	-	-

## Plastic Range

3L / 8.5L / 11.5L / 30L / 50L / 60L



Imagery not to scale

Available in small, medium and large sizes, the Clinisafe® container range offers a solution for the disposal of clinical waste within most healthcare settings.

**Key Features:** ☒ Recycled Material ☒ Hermetically Sealed ☒ Waterproof

Clinisafe® Plastic Range				Orange		Purple		Blue		Red	
Range	Size	Qty	Description	Code	NPC Code	Code	NPC Code	Code	NPC Code	Code	NPC Code
Small	3L	20		-	-	-	-	-	-	50102415	276745
Med	8.5L	20		51912410	277766	-	-	-	-	51912415	276776
Med	11.5L	20		51922410	276783	-	-	-	-	51922415	-
Large	30L	10	Flat Lid	-	-	51850420	-	51850220	-	51850415	284306
Large	30L	10	Porthole Lid	-	-	51850422	-	-	-	-	-
Large	50L	10	Flat Lid	-	-	51840420	-	-	-	51840415	-
Large	50L	10	Porthole Lid	-	-	51840422	-	-	-	-	-
Large	60L	10	Flat Lid	-	-	51860420	276806	51860220	276868	51860415	-
Large	60L	10	Porthole Lid	-	-	-	-	-	-	-	-

clinisafe®

## Cardboard Range

5L / 12L / 25L / 50L



Imagery not to scale

The Clinisafe® cardboard range offers a sustainable solution to the management of clinical waste.

**Key Features:** ✓ FSC-approved Cardboard ✓ Easy to Assemble ✓ Arrives Flat-packed for Easy Storage

Clinisafe® Cardboard Range			Orange		Purple		Blue		Gypsum	
Size	Qty	Description	Code	NHS Code	Code	NHS Code	Code	NHS Code	Code	NHS Code
5L	10		20053001	-	-	-	20055001	-	-	-
12L	10		20503001	-	-	-	20505001	-	-	-
25L	10		21003001	269426	-	-	21005001	271030	21007001	-
50L	10		21503001	-	21502001	-	21505001	271023	21507001	270415

5L, 12L, and 25L are available in pallets of 800.

50L is available in pallets of 400.



## Accessories

### Brackets / Stands / Trays / Trolleys

#### Brackets

Description	Qty	Code
0.6L Wall Bracket	10	42500000
7.5L Wall Loop	1	42750000
8L/9L/11L/13L Wall Loop	1	42150000
24L/30L Wall Loop	1	42200000
2L/3L Backboard Bracket	1	42870000
5L/7L Backboard Bracket	1	42720000
8L/9L/11L/13L Backboard Bracket	1	42410000
24L/30L Backboard Bracket	1	42210000
SafeClip Wall Bracket	10	45000000

#### Stands

Description	Qty	Code
0.6L Stand	10	41510000
1L Stand	52	42530000
1L/1.8L/2L/3L Stand (Glue Tabs)	26	42050000
5L/7L Stand	10	42570000
8L/9L/11L/13L Stand	20	42070000
High 9L Stand	10	42520004
25L Clinisafe® Cardboard Plinth	1	44100004
50L Short Clinisafe® Cardboard Plinth	1	44100001
50L Tall Clinisafe® Cardboard Plinth	1	44100005

#### Trays

Description	Qty	Code
NPSD Tray	10	43000000
Triple Tray	10	43010000
Blood Tray	12	43050000

#### Trolleys

Description	Qty	Code
Procedure Trolley	1	44000000
XL Trolley	1	44710000
Pedestal Trolley	1	44800000
SafeClip Trolley Adapter – Round 25mm	5	45000003
SafeClip Trolley Adapter – Square	5	45000004
SafeClip Trolley Adapter – Medi-Rail*	5	45000005


#### Additional

Description	Qty	Code
6L Community Specimen Bag	1	50010000
19L Community Specimen Bag	1	50020000
Absorbent Pads C135 – 150 x 150mm	2000	026100
Clinisafe® Bags	2000	020197



Imagery not to scale

\*Additional Medi-Rail adapter sizes available. Please contact your Key Account Manager for more information.



Scan me for more information about our pioneering  
Clinical Waste Management solution.

Safer today for a sustainable tomorrow.

sharpsafe<sup>®</sup> clinisafe<sup>®</sup>

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK  
+44 (0)1204 529494 [hello@vernacare.com](mailto:hello@vernacare.com) [vernacare.com](http://vernacare.com)

Vernacare  
clinical waste  
management