

CONTI[®]

CLEANSING DRY WIPES



Introducing... Conti[®] Cleansing Dry Wipes

Versatile continence care and bathing solutions
developed with personalised care in mind.

CONTI®

CLEANSING DRY WIPES



A clear simplified range

Developed and manufactured in the UK, our range provides healthcare professionals with products that support personalised care planning, whilst delivering best practice.

Our streamlined range ensure customers have full access to our most popular well recognised products, which are simpler to identify, and easier to promote.

CONTI® CLEANSING DRY WIPES LITE	CONTI® CLEANSING DRY WIPES SOFT	CONTI® CLEANSING DRY WIPES COTTON SOFT	CONTI® CLEANSING DRY WIPES SUPER SOFT	CONTI® CLEANSING DRY WIPES MACERATABLE	CONTI® CLEANSING DRY WIPES FLUSHABLE
ABSORBENCY LEVEL STRONG AND DURABLE	ABSORBENCY LEVEL STRONG AND DURABLE	ABSORBENCY LEVEL COTTON BLEND	ABSORBENCY LEVEL GENTLE TO SKIN	ABSORBENCY LEVEL CONTAINS CELLULOSE BASED FIBRES	ABSORBENCY LEVEL PLASTIC FREE WIPES
Conti® Lite combines strength and performance in a lightweight material. This everyday wipe provides a cost-effective, reliable option, for a range of care tasks.	Conti® Soft provides strength and softness to support daily cleansing needs. Soft enough for use on fragile and sensitive skin.	Conti® Cotton Soft offers a cotton and viscose enriched material for a more textile feel. Designed specifically for the effective cleansing of heavily soiled skin.	Conti® Super Soft is highly absorbent, ensuring skin is left dry without the need for vigorous rubbing and wiping. Its softness makes it ideal for cleansing all areas of the body, including fragile and sensitive skin.	Conti® Washcloth provides a more environmentally friendly alternative. Designed to disperse once disposed of in a hospital macerator, this option supports a reduction in blockage and costly callouts.	Conti® Flushable provides an alternative solution for post toileting, with the added convenience of quick, safe and hygienic disposal in either a toilet or macerator. Large and durable, this wipe is ideal for assisted toileting.

Product details

The refined, rebranded range is available to order.

Product Code	Description	Material Type	Materials Weight (GSM)	Size	Wipe Size	Wipes Per Pack	Case Qty
CLW010	Conti® Cleansing Dry Wipes Lite	100% Synthetic Thermobonded	22	Small	24 x 18cm	100	32
CLW110	Conti® Cleansing Dry Wipes Lite	100% Synthetic Thermobonded	22	Large	30 x 28cm	100	32
CSW110	Conti® Cleansing Dry Wipes Soft	Synthetic/Viscose Thermobonded	26	Large	30 x 28cm	100	32
CPV110	Conti® Cleansing Dry Wipes Cotton	Synthetic/Cotton/Viscose Thermobonded	36	Large	30 x 28cm	100	20
CXW010	Conti® Cleansing Dry Wipes Super Soft	Polyester/Viscose Spunlace	40	Small	30 x 22cm	100	16
CXW110	Conti® Cleansing Dry Wipes Super Soft	Polyester/Viscose Spunlace	40	Large	30 x 28cm	100	16
CWW075P	Conti® Cleansing Dry Wipes Maceratable	Cellulose blend Airlaid	50	Large	30 x 28cm	75	16
CF050	Conti® Cleansing Dry Wipes Flushable	Cellulose blend Wetlaid	55	Large	22 x 24cm	50	18

Vernacare



Vernacare Limited, 1 Western Avenue, Matrix Park, Chorley, Lancashire, PR7 7NB, United Kingdom.
T: +44 (0)1772 299900 F: +44 (0)1772 299901 E: info@vernacare.com www.vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
Conti® is a registered trademark of Verna Group © 2023 PH536 REV00

Scan for more info

