

International Animal Healthcare Catalogue



Vernacare
Veterinary

Contents

- 1 **Animalintex® Veterinary Dressing**
- 1 **Animalintex® Hoof Dressing**
- 2 **Gamgee® Gauze and Non-Woven Roll**
- 3 **Veterinary Cotton Wool Roll**
- 3 **Vetalintex® Hydrogel Dressing**
- 4 **Skintact Low Adherent Dressing**
- 5 **Equiwrap® Cohesive Bandage**
- 6 **Flexoplast® Elastic Adhesive Bandage**
- 6 **Fast Aid Orthopaedic Padding Bandage**
- 7 **Stayform® Dressing Retention Bandage**
- 7 **Fast Aid® Zinc Oxide Tape**



Animalintex® Veterinary Dressing

The convenient multi-layered veterinary dressing that can be used as a hot wet poultice, cold wet poultice or as a dry dressing. Animalintex® is ideal for a wide variety of applications and a key staple for any first aid kit.

Key Features & Benefits

- Certified quality of manufacture
- Trusted, safe and stringently tested
- Manufactured in the UK
- Tried and tested through decades of use



| Product Code | Description | Size | Unit of Sale |
|--------------|---------------------------------------|---------------|------------------|
| 5704 | Animalintex Dressing – EU Only | 209mm x 406mm | 6 x 10 dressings |
| 8703 | Animalintex Dressing – Export | 209mm x 406mm | 6 x 10 dressings |
| 8704 | Animalintex Dressing – Australia Only | 209mm x 406mm | 6 x 10 dressings |

Animalintex® Hoof Dressing

The convenient multi-layered veterinary dressing that can be used as a hot wet poultice, cold wet poultice or as a dry dressing. Animalintex® is ideal for a wide variety of applications and a key staple for any first aid kit.

Animalintex® hoof dressing is pre-cut to shape for economy and convenience.

Key Features & Benefits

- Certified quality of manufacture
- Trusted, safe and stringently tested
- Manufactured in the UK
- Tried and tested through decades of use



| Product Code | Description | Size | Unit of Sale |
|--------------|---------------------------|---------|------------------|
| 8719 | Animalintex Hoof – Export | 140mm Ø | 12 x 3 dressings |



Gamgee® Gauze and Non-Woven Roll

The original Gamgee® from Robinson Animal Healthcare consists of a thick layer of high quality, absorbent cotton enclosed in either traditional gauze or a non-woven cover. Insulates and cushions the wound, protecting it from external trauma.

Key Features & Benefits

- Supplied in 500g rolls
- Available in 30cm and 45cm widths
- Option of gauze or non-woven cover
- Highly absorbent secondary wound dressing
- Provides insulation, padding and protection to minimise wound trauma
- Ideal for general swabbing
- Sterile Gamgee® cut pieces are also available



| Product Code | Description | Width | Unit of Sale |
|--------------|---------------------------------------|-------|--------------|
| 6556 | Veterinary Gamgee Non-Woven Roll 500g | 30cm | 30 rolls |
| 6557 | Veterinary Gamgee Non-Woven Roll 500g | 45cm | 24 rolls |
| 8709 | Veterinary Gamgee Gauze Roll 500g | 30cm | 30 rolls |
| 8710 | Veterinary Gamgee Gauze Roll 500g | 45cm | 24 rolls |

Veterinary Cotton Wool Roll

Premium quality cotton wool rolls, manufactured from 100% pure cotton that is soft and highly absorbent. Suitable as a cushioning and highly absorbent secondary wound layer, providing protection and insulation or is also ideal for topical cleansing and swabbing.

Key Features & Benefits

- 100% natural and sustainable
- Soft and highly absorbent
- Suitable as a cushioning and absorbent secondary wound layer
- Ideal for topical cleansing, even around delicate areas
- Available in either a 350g or 500g roll
- Supplied in a convenient drawstring bag



| Product Code | Description | Weight | Unit of Sale |
|--------------|-----------------------------|--------|--------------|
| 6555 | Veterinary Cotton Wool Roll | 350g | 12 rolls |
| 6560 | Veterinary Cotton Wool Roll | 500g | 12 rolls |

Vetalintex® Hydrogel Dressing

A clear, viscose starch based hydrogel dressing. The hydrogel is supplied sterile and has a gentle cleansing action.

Key Features & Benefits

- Promotes wound healing by maintaining optimum hydration for moist wound cleansing
- Excellent, gentle cleansing action
- Secure screw cap with soft tube for ease of application
- Presented sterile for optimum infection control



| Product Code | Description | Weight | Unit of Sale |
|--------------|---------------------------------|----------|---------------|
| 6563 | Vetalintex Sterile Hydrogel 15g | 15g tube | 20 x 15 tubes |

Skintact Low Adherent Dressing

A high quality, sterile, perforated film dressing that is used as a low adherent wound contact layer. The dressing is the same on both sides, preventing the incorrect side being used on the wound. It consists of two layers of perforated film with an inner absorbent layer. The perforations allow drainage of exudates to the absorbent layer.

Key Features & Benefits

- Low adherent wound contact layer to aid removal and minimise trauma to the wound
- Double sided, can be used either side
- Pads, cushions and protects wounds
- Supplied in individual sterile packs and dispensing boxes



| Product Code | Description | Size | Unit of Sale |
|--------------|--------------------------------|-------------|-------------------|
| 7900 | Skintact Low Adherent Dressing | 5cm x 5cm | 6 x 100 dressings |
| 7901 | Skintact Low Adherent Dressing | 10cm x 10cm | 6 x 100 dressings |
| 7902 | Skintact Low Adherent Dressing | 10cm x 20cm | 6 x 100 dressings |





Equiwrap® Cohesive Bandage

A flexible, lightweight and highly effective cohesive bandage that is ideal for support, protection and dressing retention. Particularly suited for dressing awkward contours and areas, each bandage tears easily and accurately by hand for quick application.

Key Features & Benefits

- Flexible and lightweight
- Ideal for dressing retention and support
- Quick and easy to apply
- Sticks to itself, but not skin or hair
- Available in a wide range of colours



| Product Code | Description | Size | Unit of Sale |
|--------------|---|-------------|-------------------|
| 8660 | Equiwrap Cohesive Bandage – Black | 10cm x 4.5m | 12 bandages rolls |
| 8661 | Equiwrap Cohesive Bandage – Blue | 10cm x 4.5m | 12 bandages rolls |
| 8663 | Equiwrap Cohesive Bandage – Green | 10cm x 4.5m | 12 bandages rolls |
| 8664 | Equiwrap Cohesive Bandage – Red | 10cm x 4.5m | 12 bandages rolls |
| 8665 | Equiwrap Cohesive Bandage – White | 10cm x 4.5m | 12 bandages rolls |
| 8666 | Equiwrap Cohesive Bandage – Yellow | 10cm x 4.5m | 12 bandages rolls |
| 8668 | Equiwrap Cohesive Bandage – Purple | 10cm x 4.5m | 12 bandages rolls |
| 8670 | Equiwrap Cohesive Bandage – Neon Green | 10cm x 4.5m | 12 bandages rolls |
| 8671 | Equiwrap Cohesive Bandage – Neon Orange | 10cm x 4.5m | 12 bandages rolls |
| 8672 | Equiwrap Cohesive Bandage – Neon Pink | 10cm x 4.5m | 12 bandages rolls |
| 8667 | Equiwrap Cohesive Bandage – Mixed Colours | 10cm x 4.5m | 24 bandages rolls |

Flexoplast® Elastic Adhesive Bandage

An elastic adhesive bandage supplied with an adhesive outer face to make application easier, particularly on smaller animals. Manufactured from strong, woven elasticated fabric, each bandage has a centre line to aid application. Suitable for immobilisation, support and dressing retention.

Key Features & Benefits

- Strong, adhesive bandage
- Adhesive coated outer face and centre guideline for ease of application
- Suitable for immobilisation, support and dressing retention
- Ideal for awkward to dress areas



| Product Code | Description | Size | Unit of Sale |
|--------------|-------------------------------------|--------------|--------------|
| 7880 | Flexoplast Elastic Adhesive Bandage | 5cm x 4.5m | 12 bandages |
| 7881 | Flexoplast Elastic Adhesive Bandage | 7.5cm x 4.5m | 12 bandages |
| 7882 | Flexoplast Elastic Adhesive Bandage | 10cm x 4.5m | 12 bandages |

Fast Aid Orthopaedic Padding Bandage

A soft and protective bandage that is ideal for the cushioning and protection of bony prominences. It may be applied directly to the skin for padding splints and under adhesive bandages. Manufactured from non absorbent polyester that allows moisture to drain away from the skin whilst retaining its shape and loft, even when wet.

Key Features & Benefits

- Soft, smooth and highly conformable
- Protects bony prominences
- Ideal as under cast padding
- Easy to apply and tear
- Retains shape and loft, even when wet
- Individually packed



| Product Code | Description | Size | Unit of Sale |
|--------------|--------------------------------------|--------------|--------------|
| 7670 | Fast Aid Orthopaedic Padding Bandage | 5cm x 2.7m | 12 bandages |
| 7671 | Fast Aid Orthopaedic Padding Bandage | 7.5cm x 2.7m | 12 bandages |
| 7672 | Fast Aid Orthopaedic Padding Bandage | 10cm x 2.7m | 12 bandages |

Stayform® Dressing Retention Bandage

A high quality woven polyamide and cellulose contour bandage, designed to retain dressings securely to the wound without inhibiting movement or circulation. Thin, lightweight and comfortable in use, each bandage has non-fraying edges and is kind and gentle to the most sensitive of skin.

Key Features & Benefits

- Retains dressings securely
- Soft and conformable
- Will not restrict movement or circulation
- Thin, lightweight and comfortable in use
- Gentle and kind to delicate skin
- Supplied individually wrapped



| Product Code | Description | Size | Unit of Sale |
|--------------|-------------------------------------|------------|-----------------|
| 5680 | Stayform Dressing Retention Bandage | 5cm x 4m | 12 x 6 bandages |
| 5681 | Stayform Dressing Retention Bandage | 7.5cm x 4m | 12 x 6 bandages |
| 5682 | Stayform Dressing Retention Bandage | 10cm x 4m | 12 x 6 bandages |
| 5683 | Stayform Dressing Retention Bandage | 15cm x 4m | 6 x 6 bandages |

Fast Aid® Zinc Oxide Tape

Can be applied directly to the skin and helps to support joints by immobilising or restricting movement. The strong, tan coloured tape can also be used to secure dressings in place or to protect soft tissue and help prevent blisters.

Key Features & Benefits

- Restricts joint movement to help prevent against injury
- Helps to protect the skin from blisters
- Suitable to secure dressings in place
- Supplied in individual boxes



| Product Code | Description | Size | Unit of Sale |
|--------------|--------------------------|-------------|--------------|
| 4482 | Fast Aid Zinc Oxide Tape | 1.25cm x 5m | 12 rolls |
| 4483 | Fast Aid Zinc Oxide Tape | 2.5cm x 5m | 12 rolls |
| 4484 | Fast Aid Zinc Oxide Tape | 5cm x 5m | 12 rolls |



**Scan me to find out more about
Vernacare's products and solutions.**

Robinson Animal Healthcare Ltd
Lawn Road
Carlton-in-Lindrick
Worksop
S81 9LB, UK
Tel: +44 (0)1909 735000
Fax: +44 (0)1909 735013
www.robinsonanimalhealthcare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
© Copyright Vernacare International Limited 2025. All rights reserved. SU 1299 REV00

The logo for Vernacare Veterinary, featuring a white wave-like line above the text "Vernacare" and "Veterinary" stacked vertically.