

UK Animal Healthcare Catalogue



Contents

- 1 **Animalintex® Licensed Poultice Dressing**
- 1 **Animalintex® Licensed Hoof Poultice Dressing**
- 2 **Gamgee® Gauze and Non-Woven Roll**
- 3 **Veterinary Cotton Wool Roll**
- 3 **Vetalintex® Hydrogel Dressing**
- 4 **Skintact Low Adherent Dressing**
- 5 **Equiwrap® Cohesive Bandage**
- 6 **Flexoplast® Elastic Adhesive Bandage**
- 6 **Fast Aid Orthopaedic Padding Bandage**
- 7 **Stayform® Dressing Retention Bandage**
- 7 **Fast Aid® Zinc Oxide Tape**



Animalintex® Licensed Poultice Dressing

Animalintex® dressings are the convenient, multi-layered, absorbent veterinary poultice licensed in the UK by the Veterinary Medicines Directorate (VMD). It can be used as a hot wet poultice, cold wet poultice or dry dressing and is suitable for the treatment of infected wounds, open wounds, foot punctures, corns, bruised soles, thrush, seedy toe, abscesses, laminitis, sandcracks and thorns.

Key Features & Benefits

- Safe, licensed veterinary medicinal product
- VMD audited manufacturing facility for consistent quality and effectiveness
- Manufactured in the UK
- Can be used as a hot or cold poultice or as a dry dressing
- Tried and tested through decades of use



Marketing Authorisation No. Vm19434/4001

Product Code	Description	Size	Unit of Sale
76010	Animalintex Licensed Poultice Dressing*	209mm x 406mm	6 x 10 dressings

*Animalintex Licensed Poultice Dressing is strictly for sale, distribution and use in the UK only.

Animalintex® Licensed Hoof Poultice Dressing

Animalintex® dressings are the convenient, multi-layered, absorbent veterinary poultice licensed in the UK by the Veterinary Medicines Directorate (VMD). It can be used as a hot wet poultice, cold wet poultice or dry dressing and is ideal for the treatment of foot conditions such as punctures, corns, bruised soles, thrush, seedy toe, abscesses, laminitis, sandcracks and thorns.

Animalintex® hoof dressing is pre-cut to shape for economy and convenience.

Key Features & Benefits

- Safe, licensed veterinary medicinal product
- VMD audited manufacturing facility for consistent quality and effectiveness
- Manufactured in the UK
- Can be used as a hot or cold poultice or as a dry dressing
- Tried and tested through decades of use



Marketing Authorisation No. Vm19434/4000

Product Code	Description	Size	Unit of Sale
8717	Animalintex Licensed Hoof Poultice Dressing*	140mm Ø	12 x 3 dressings

*Animalintex Licensed Hoof Poultice Dressing is strictly for sale, distribution and use in the UK only.



Gamgee® Gauze and Non-Woven Roll

The original Gamgee® from Robinson Animal Healthcare consists of a thick layer of high quality, absorbent cotton enclosed in either traditional gauze or a non-woven cover. Insulates and cushions the wound, protecting it from external trauma.

Key Features & Benefits

- Supplied in 500g rolls
- Available in 30cm and 45cm widths
- Option of gauze or non-woven cover
- Highly absorbent secondary wound dressing
- Provides insulation, padding and protection to minimise wound trauma
- Ideal for general swabbing
- Sterile Gamgee® cut pieces are also available



Product Code	Description	Width	Unit of Sale
6556	Veterinary Gamgee Non-Woven Roll 500g	30cm	30 rolls
6557	Veterinary Gamgee Non-Woven Roll 500g	45cm	24 rolls
8709	Veterinary Gamgee Gauze Roll 500g	30cm	30 rolls
8710	Veterinary Gamgee Gauze Roll 500g	45cm	24 rolls

Veterinary Cotton Wool Roll

Premium quality cotton wool rolls, manufactured from 100% pure cotton that is soft and highly absorbent. Suitable as a cushioning and highly absorbent secondary wound layer, providing protection and insulation or is also ideal for topical cleansing and swabbing.

Key Features & Benefits

- 100% natural and sustainable
- Soft and highly absorbent
- Suitable as a cushioning and absorbent secondary wound layer
- Ideal for topical cleansing, even around delicate areas
- Available in either a 350g or 500g roll
- Supplied in a convenient drawstring bag



Product Code	Description	Weight	Unit of Sale
6555	Veterinary Cotton Wool Roll	350g	12 rolls
6560	Veterinary Cotton Wool Roll	500g	12 rolls

Vetalintex® Hydrogel Dressing

A clear, viscose starch based hydrogel dressing. The hydrogel is supplied sterile and has a gentle cleansing action.

Key Features & Benefits

- Promotes wound healing by maintaining optimum hydration for moist wound cleansing
- Excellent, gentle cleansing action
- Secure screw cap with soft tube for ease of application
- Presented sterile for optimum infection control



Product Code	Description	Weight	Unit of Sale
6563	Vetalintex Sterile Hydrogel 15g	15g tube	20 x 15 tubes

Skintact Low Adherent Dressing

A high quality, sterile, perforated film dressing that is used as a low adherent wound contact layer. The dressing is the same on both sides, preventing the incorrect side being used on the wound. It consists of two layers of perforated film with an inner absorbent layer. The perforations allow drainage of exudates to the absorbent layer.

Key Features & Benefits

- Low adherent wound contact layer to aid removal and minimise trauma to the wound
- Double sided, can be used either side
- Pads, cushions and protects wounds
- Supplied in individual sterile packs and dispensing boxes



Product Code	Description	Size	Unit of Sale
7900	Skintact Low Adherent Dressing	5cm x 5cm	6 x 100 dressings
7901	Skintact Low Adherent Dressing	10cm x 10cm	6 x 100 dressings
7902	Skintact Low Adherent Dressing	10cm x 20cm	6 x 100 dressings





Equiwrap® Cohesive Bandage

A flexible, lightweight and highly effective cohesive bandage that is ideal for support, protection and dressing retention. Particularly suited for dressing awkward contours and areas, each bandage tears easily and accurately by hand for quick application.

Key Features & Benefits

- Flexible and lightweight
- Ideal for dressing retention and support
- Quick and easy to apply
- Sticks to itself, but not skin or hair
- Available in a wide range of colours



Product Code	Description	Size	Unit of Sale
8660	Equiwrap Cohesive Bandage – Black	10cm x 4.5m	12 bandages rolls
8661	Equiwrap Cohesive Bandage – Blue	10cm x 4.5m	12 bandages rolls
8663	Equiwrap Cohesive Bandage – Green	10cm x 4.5m	12 bandages rolls
8664	Equiwrap Cohesive Bandage – Red	10cm x 4.5m	12 bandages rolls
8665	Equiwrap Cohesive Bandage – White	10cm x 4.5m	12 bandages rolls
8666	Equiwrap Cohesive Bandage – Yellow	10cm x 4.5m	12 bandages rolls
8668	Equiwrap Cohesive Bandage – Purple	10cm x 4.5m	12 bandages rolls
8670	Equiwrap Cohesive Bandage – Neon Green	10cm x 4.5m	12 bandages rolls
8671	Equiwrap Cohesive Bandage – Neon Orange	10cm x 4.5m	12 bandages rolls
8672	Equiwrap Cohesive Bandage – Neon Pink	10cm x 4.5m	12 bandages rolls
8667	Equiwrap Cohesive Bandage – Mixed Colours	10cm x 4.5m	24 bandages rolls

Flexoplast® Elastic Adhesive Bandage

An elastic adhesive bandage supplied with an adhesive outer face to make application easier, particularly on smaller animals. Manufactured from strong, woven elasticated fabric, each bandage has a centre line to aid application. Suitable for immobilisation, support and dressing retention.

Key Features & Benefits

- Strong, adhesive bandage
- Adhesive coated outer face and centre guideline for ease of application
- Suitable for immobilisation, support and dressing retention
- Ideal for awkward to dress areas



Product Code	Description	Size	Unit of Sale
7880	Flexoplast Elastic Adhesive Bandage	5cm x 4.5m	12 bandages
7881	Flexoplast Elastic Adhesive Bandage	7.5cm x 4.5m	12 bandages
7882	Flexoplast Elastic Adhesive Bandage	10cm x 4.5m	12 bandages

Fast Aid Orthopaedic Padding Bandage

A soft and protective bandage that is ideal for the cushioning and protection of bony prominences. It may be applied directly to the skin for padding splints and under adhesive bandages. Manufactured from non absorbent polyester that allows moisture to drain away from the skin whilst retaining its shape and loft, even when wet.

Key Features & Benefits

- Soft, smooth and highly conformable
- Protects bony prominences
- Ideal as under cast padding
- Easy to apply and tear
- Retains shape and loft, even when wet
- Individually packed



Product Code	Description	Size	Unit of Sale
7670	Fast Aid Orthopaedic Padding Bandage	5cm x 2.7m	12 bandages
7671	Fast Aid Orthopaedic Padding Bandage	7.5cm x 2.7m	12 bandages
7672	Fast Aid Orthopaedic Padding Bandage	10cm x 2.7m	12 bandages

Stayform® Dressing Retention Bandage

A high quality woven polyamide and cellulose contour bandage, designed to retain dressings securely to the wound without inhibiting movement or circulation. Thin, lightweight and comfortable in use, each bandage has non-fraying edges and is kind and gentle to the most sensitive of skin.

Key Features & Benefits

- Retains dressings securely
- Soft and conformable
- Will not restrict movement or circulation
- Thin, lightweight and comfortable in use
- Gentle and kind to delicate skin
- Supplied individually wrapped



Product Code	Description	Size	Unit of Sale
5680	Stayform Dressing Retention Bandage	5cm x 4m	12 x 6 bandages
5681	Stayform Dressing Retention Bandage	7.5cm x 4m	12 x 6 bandages
5682	Stayform Dressing Retention Bandage	10cm x 4m	12 x 6 bandages
5683	Stayform Dressing Retention Bandage	15cm x 4m	6 x 6 bandages

Fast Aid® Zinc Oxide Tape

Can be applied directly to the skin and helps to support joints by immobilising or restricting movement. The strong, tan coloured tape can also be used to secure dressings in place or to protect soft tissue and help prevent blisters.

Key Features & Benefits

- Restricts joint movement to help prevent against injury
- Helps to protect the skin from blisters
- Suitable to secure dressings in place
- Supplied in individual boxes



Product Code	Description	Size	Unit of Sale
4482	Fast Aid Zinc Oxide Tape	1.25cm x 5m	12 rolls
4483	Fast Aid Zinc Oxide Tape	2.5cm x 5m	12 rolls
4484	Fast Aid Zinc Oxide Tape	5cm x 5m	12 rolls



**Scan me to find out more about
Vernacare's products and solutions.**

Robinson Animal Healthcare Ltd
Lawn Road
Carlton-in-Lindrick
Worksop
S81 9LB, UK
Tel: +44 (0)1909 735000
Fax: +44 (0)1909 735013
www.robinsonanimalhealthcare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.
© Copyright Vernacare International Limited 2025. All rights reserved. SU1300 REV00

The logo for Vernacare Veterinary, featuring a white wavy line above the word "Vernacare" in a bold, sans-serif font, with the word "Veterinary" in a smaller, sans-serif font below it.